

## Lab Protocols

### Lab Protocols

- Bisulfite treatment protocol [1]
- Updated method for the isolation of pancreatic islets [2]
- Immunohisto Slide Processing [3]
- Chromatin IP General Protocol [4]
- A rapid micro chromatin immunoprecipitation assay (iChIP) [5]

Image not found

[islet](https://bhushanlab.med.ucla.edu/images/islet_sm.jpg)

[Contact Us](#)  
[UCSF Main Site](#)

© 2015 The Regents of the University of California

---

**Source URL:** <https://bhushanlab.ucsf.edu/protocols>

#### Links

[1] [http://bhushanlab.med.ucla.edu/downloads/Bisulfite\\_treatment\\_protocol.pdf](http://bhushanlab.med.ucla.edu/downloads/Bisulfite_treatment_protocol.pdf)

[2] [http://bhushanlab.med.ucla.edu/downloads/Georgia\\_Islet\\_Isolation.pdf](http://bhushanlab.med.ucla.edu/downloads/Georgia_Islet_Isolation.pdf)

[3] [http://bhushanlab.med.ucla.edu/downloads/Immunohisto\\_slide\\_processing.pdf](http://bhushanlab.med.ucla.edu/downloads/Immunohisto_slide_processing.pdf)

[4] [http://bhushanlab.med.ucla.edu/downloads/ChIP\\_Protocol.pdf](http://bhushanlab.med.ucla.edu/downloads/ChIP_Protocol.pdf)

[5] [http://bhushanlab.med.ucla.edu/downloads/Micro\\_ChIP\\_Nature\\_Protocols.pdf](http://bhushanlab.med.ucla.edu/downloads/Micro_ChIP_Nature_Protocols.pdf)